

Sirius Optics Unit 1 26 Darnick Street Underwood, Qld 4119 Opening Hours

10am-5:30pm Mon-Fri 9am-2pm Sat Phone: 07 3423 2355 www.sirius-optics.com.au

# saxon Microscope Immersion Oil

## AUD \$19.95

#### **Product Images**



### **Short Description**

For use in light microscopy.

In light microscopy, oil immersion is a technique used to increase the resolving power of a microscope. This is achieved by immersing both the objective lens and the specimen in a transparent oil of high refractive index, thereby increasing the numerical aperture of the objective lens. NOTE: Immersion Oil is supplied in 5ml qty.

#### **Additional Information**

Specifications

No